





# Kosmos Chemical Compatibility

The cleaning and disinfection agents listed below have been evaluated for compatibility with Kosmos device materials. Refer to the agent manufacturer instructions for cleaning and disinfecting procedures and contact conditions.

Cleaning and/or Disinfectant Agent	Active Ingredient	Kosmos Bridge	Kosmos Probes
Warm soapy water	None	✓	✓
PDI Healthcare Sani-Cloth® Plus Germicidal Disposable Wipes	n-Alkyl dimethyl ethylbenzyl ammonium chlorides (0.125%), n-Alkyl dimethyl Benzyl Ammonium Chlorides (0.125%), Isopropanol (10-20%), 2-Butoxyethanol (<5%)	✓	✓
¹PDI Healthcare Super Sani-Cloth® Germicidal Disposable Wipes	n-Alkyl dimethyl ethylbenzyl ammonium chlorides (0.25%), n-Alkyl dimethyl Benzyl Ammonium Chlorides (0.25%), Isopropyl Alcohol (Isopropanol) (55.5%)	✓	✓
¹PDI Healthcare Sani-Cloth® Bleach Germicidal Disposable Wipes	Sodium Hypochlorite (0.63%)	✓	✓
PDI Healthcare Sani-Cloth™ Prime Wipes	Didecyl Dimethyl Ammonium Chloride (0.61%), Ethyl Alcohol (27.30%), Isopropyl Alcohol (28.70%), Other ingredients (43.39%)	N/A	✓
Sani-HP1™ Germicidal Wipes	Hydrogen Peroxide (0.58%), Other ingredients (99.42%)	N/A	✓
³Metrex CaviWipes™ (KavoWipes)	Diisobutylphenoxyethoxyethyl dimethylbenzyl ammonium chloride (0.28%), Isopropanol (17.2%)	✓	✓
¹Metrex CaviCide™	Diisobutylphenoxyethoxyethyl dimethylbenzyl ammonium chloride (0.28%), Isopropanol (17.2%)	✓	✓
¹Metrex CaviCide1™	Didecyl dimethyl ammonium Chloride (0.76%), Ethanol (7.50%), Isopropanol (15%)	✓	✓
ASP Cidex® OPA	ortho-Phthalaldehyde (0.55%), Dipotassium Hydrogen phosphate, Potassium Dihydrogen Phosphate, Benzotriazole, citric acid, D&C green dye #5, N-(Hydroxyethyl)-ethylenediaminetriacetic acid (HEDTA) (99.45%)	N/A	✓
ASP Enzo® Enzymatic Detergent	Subtilisins (<5%)	N/A	✓
Hydrogen Peroxide 3% solution	Hydrogen Peroxide (3%), water (97%)	✓	✓
Kyorin Pharmaceutical Milton	Sodium Hypochlorite (1.1%)	✓	✓
Tristel Duo ULT™	Sodium Chlorite (0.5%), Citric Acid (5.26%), Chlorine Dioxide (0.032%)	N/A	✓
Trophon2 and Trophon3	Hydrogen Peroxide (35%)	N/A	✓

Key	
<sup>1</sup> : Product on List	Disinfectants for Use Agent SARS-CoV-2
<sup>2</sup> : Per Metrex released note	Approved for Covid 19
	Compatible
	Compatibility data not available



Scan the QR code to access  
additional product resources  
and support materials

## Customer Support



**EchoNous, Inc.**  
8310 154th Ave. NE, Building B, Suite 200  
Redmond, WA 98052

Call 1-844-854-0800  
Email [support@echonous.com](mailto:support@echonous.com)

Refer to user guide before use.

⚠ **Caution:** Federal (United States) law restricts this device to sale by or on the order of a physician.