

Summary


This technical bulletin is for installing the EchoNous Bladder Probe replacement cable. This EchoNous cable kit part number is **P006376 or R006376 (Refurbished)**.

1. Remove the existing cable from the EchoNous Bladder Probe by unscrewing (counterclockwise) the gold ring at the base of the probe handle by hand and unplugging the cable from the probe body.

DO NOT USE PLIERS OR VICE GRIPS AS THIS WILL DAMAGE THE PROBE.

If you are not able to remove the cable by hand, please contact EchoNous Customer Service at 844-854-0800.

2. Remove the replacement cable from the bag and verify that the serial number on the blue tag matches the one you are replacing.
3. Before threading the cable's locking nut onto the probe, make sure the o-ring at the end of the cable is seated against USB connector **(A)**, and that the black probe body gasket is in place **(B)**.
4. Insert the new cable into the USB-C connector of the probe **(C)**, and thread the gold locking nut onto the bladder probe and **finger tighten only (D)**.
5. Plug the EchoNous Bladder Probe into the hub or directly into the tablet and the test the scanner. Select Bladder Volume from the home screen, select a patient type and run a scan. This step only checks that there is communication between the tablet and the probe. It does not check the accuracy of the probe.



If you have any questions regarding the installation process, contact EchoNous Customer Support at 844-854-0800.