

Quick Start Guide

Kosmos and 19Labs

Inviting a 2nd Camera to your Mentorship Sessions

How to give the instructor a "birds-eye view" of your patient/hand/probe using a 2nd camera. Please complete the following steps during the connection phase of you of your sessions if possible.




Figure 1.1




Figure 1.2

1. Ensure you have a Laptop with a camera (built-in or USB webcam) or Smartphone with a camera and Zoom installed with an email address associated with the device, connected to WiFi, and the ability to position the camera to focus on your patient.
2. Set up your computer, webcam, or phone to face the patient area for scanning. For laptops with flip-lid screens, position them on a shelf and tilt the camera/screen downward toward the patient. Be creative with other setups as needed.
3. Start your mentorship connection on your Kosmos tablet in the usual manner following steps 1-3 on the ACP Mentorship Quick Start Guide.
4. After connecting to your mentor's clinic through the ACP Mentorship tablet, locate the "INVITE" icon in the top toolbar alongside the camera icons. Click on it to send the Zoom invite via email to your second device with the camera. Alternatively, if you're using your phone as the camera, utilize the SMS option to send the invite. (See figure 1.1)
5. Once your second camera device receives the Zoom link, simply click to join the Zoom session as a new member. This allows your mentor to view both cameras: the laptop/ phone and the Kosmos tablet's camera.
6. If you have a USB webcam connected to your laptop along with its built-in camera, you may need to switch between camera options. Simply click on the "up carrot" on the Video button (Figure 1.2) to do so.
7. Now complete the rest of the steps (Step 4) on the ACP Quick Start Guide to connect to your Mentorship Session.

Contact

Customer Support

Call 1-844-854-0800

Email service@echonous.com



D011591_Rev A



19Labs is a third party provided software by EchoNous. EchoNous therefore disclaims any and all liability, including any express or implied warranties, whether oral or written, for 19Labs services. By purchasing 19Labs' products or services, customers acknowledge that no representation has been made by EchoNous as to the fitness of the third-party services for the intended purpose. All 19Labs software purchased by customer is governed by the Terms of Use, warranties and disclaimers published by 19Labs. Any questions regarding the 19Labs Terms of Use should be brought directly to 19Labs at: info@19labs.com

©2024 EchoNous, Inc. EchoNous and KOSMOS are registered trademarks of EchoNous, Inc.