

# The EchoNous Carrying Case

Hand-carried stand and  
storage for the EchoNous  
System



## Features



### Dual Purpose

Easily connect a probe of the EchoNous System to the tablet for live scanning.



### Rigid Structure

Built-in stiffeners and the internal divider panel provide stability so that the EchoNous carrying case won't lean or collapse.



### Convertible front cover

With just a simple flip, the protective front cover can fold underneath the carrying case to provide a stable base and reveal the tablet that is securely attached for viewing while scanning.



### Ample Storage

Open the back panel to reveal internal storage for supplies and mesh pockets to hold the probe, charger and e-box.



### Adaptive Design

The tablet frame can be adjusted to accommodate different sizes for future upgradeability.

## Specifications

### Weight

- 800 grams (excluding tablet, e-box and probes)

### Dimensions

- 7.67" D x 12.5" W x 8.8" H

### Protective Materials

- Non-skid base materials
- Rubber front cover
- Internal stiffeners built into main compartment

### Storage

- Exterior mesh pocket
- Interior storage compartment
- Interior mesh pockets