

EchoNous AI Station

For EchoNous Bladder and Vein



FEATURES



Versatile and Agile

Easily roll the AI Station from the front or rear handle and adjust the tablet mount 180° for an optimal view of the screen



Rugged Tablet Mount

Designed to be durable and withstand fast paced environments by including built-in USB cable management and a sturdy mount hinge



Intuitive Design

With just a push of a quick release button, the rear storage bin can be removed for cleaning

CONTACT

www.EchoNous.com

8310 154th Avenue NE Bldg B, Suite 200
Redmond, WA 98052 United States
Tel 844-854-0800

SPECIFICATIONS

Weight	• 24.5 kg
Dimensions	• 137 H x 41 W x 59 D cm
Base Dimensions	• 36 W x 46 D cm
Storage	• Quick release, rear mounted storage bin • Holder for ultrasound gel
Cable Management	• Central track with covered panel
Casters	• 5" locking casters
Component Adjustability	• Handle bar and tablet mount can be installed at preferred height
Mounting Bracket	• Swivel left and right (+/- 90°) • Tilt up and down (+30°/-15°)
Durability	• Heavy duty base, column and hinge • Protection for USB connectors
Warranty	• 5-year warranty