
Product: Kosmos Tablet G1

Date: 8/10/2026

We're pleased to announce the release of Kosmos Tablet G1, a rugged Android tablet for the Kosmos ultrasound system. Custom-engineered through an exclusive partnership with Getac, a global leader in ruggedized mobile technology, Kosmos Tablet G1 is a variant of the Getac ZX10 tablet built for and only available through EchoNous.

This major update creates a seamless, out-of-the-box experience, allowing the user to boot up and start scanning in the most demanding environments.

Kosmos Tablet G1

- Built with Getac from a variant of the ZX10 tablet
- The lightest 10" fully-rugged Android tablet on the market
- Certified MIL-STD 810H with a 6-foot drop rating
- IP66-rated against dust and high-pressure water
- Operable from -20F to 145F
- Sunlight-readable Lumibond display
- Wi-Fi 6E
- DICOM TLS Capable
- Swappable Batteries for extended scanning

System Overview

The Kosmos Tablet G1 includes a wide range of features and capabilities allowing customers to create a customized device based on their specific needs. Please note, not all features and products are available in all regions. Please speak with your local EchoNous representative to check availability.

Compatible Hardware/Transducers

- Kosmos Torso-One USB*
- Kosmos Link
- Kosmos Tablet G1 Handstrap/Kickstand
- Kosmos Tablet G1 Dual Battery Charger
- Kosmos Lexsa
- Kosmos Tablet G1 High Capacity Battery
- Kosmos Tablet G1 Standard Battery

Available Presets*

- Kosmos Bladder
- OB & GYN
- TCD
- Torso-One Core Presets (Heart, Lung, Abdomen)
- Lexsa Core Presets (MSK, Nerve, Vascular, Lung)

Capabilities

- B-Mode/2D
- Color Doppler
- Pulse Wave Doppler
- Tissue Doppler Imaging
- M-Mode
- Color Power Doppler (Lexsa only)
- Continuous-Wave Doppler (Torso-One only)

Automated Workflows

- AI FAST
- Auto Doppler
- Auto Preset
- Auto EF
- Auto VTI Trace

If you have any questions regarding feature activation, contact EchoNous Customer Service by telephone at (844) 854-0800 or via email at service@echonous.com.

*Important Notifications

- This list includes purchasable features. Please note, not all features and products are available in all regions. Please speak with your local EchoNous representative to check availability.
- Kosmos Torso-One B (for Bridge devices) is not compatible with Kosmos Tablet G1
- Kosmos EDU are not currently available for the Kosmos Tablet G1