

Purchasable Software Features for Kosmos on Bridge, Android, and iOS

Date: August 23rd, 2023

Distribution: Internal and external (not for general distribution).

This document outlines the available features when a Kosmos software version is purchased. Please contact EchoNous support for any questions regarding software availability.

Kosmos Bridge: Purchasable Software Features

Software Version	Features
v3.0	Torso/Torso-One Features: <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • Trio • AI-assisted EF Workflow • PW Doppler • CW Doppler
v3.2	Torso/Torso-One Features: <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • Trio • AI-assisted EF Workflow • PW Doppler • CW Doppler
v4.1	Torso/Torso-One Features: <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • Trio • AI-assisted EF Workflow • PW Doppler • CW Doppler

	Lexsa Features: <ul style="list-style-type: none"> • B-mode • M-mode
v4.1.3 (Canada)	Torso/Torso-One Features: <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • AI-assisted EF Workflow • PW Doppler • CW Doppler Lexsa Features: <ul style="list-style-type: none"> • B-mode • M-mode
v4.3	Integrated Software Features: <ul style="list-style-type: none"> • Kosmos Hub compatibility and support • Screen mirroring (Miracast) function • VPN configurations Torso/Torso-One Features: <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • Trio • AI-assisted EF Workflow • PW Doppler • CW Doppler Lexsa Features: <ul style="list-style-type: none"> • B-mode • M-mode
v4.5	Integrated Software Features: <ul style="list-style-type: none"> • FIPS 140-2 Compliance (US only) Torso/Torso-One Features: <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • Trio • AI-assisted EF Workflow

	<ul style="list-style-type: none"> • PW Doppler • CW Doppler Lexsa Features: <ul style="list-style-type: none"> • B-mode • M-mode
v4.6	Integrated Software Features: <ul style="list-style-type: none"> • Kosmos UP (US only) Torso/Torso-One Features: <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • Trio • AI-assisted EF Workflow • PW Doppler • CW Doppler Lexsa Features: <ul style="list-style-type: none"> • B-mode • M-mode
v6.0	Torso/Torso-One Features: <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • Trio • AI-assisted EF Workflow • PW Doppler • CW Doppler Lexsa Features: <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler
v7.0	Integrated Software Features: <ul style="list-style-type: none"> • Us2.ai compatibility Torso/Torso-One Features: <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • Trio

	<ul style="list-style-type: none"> • AI-assisted EF Workflow • PW Doppler • CW Doppler • TDI <p>Lexsa Features:</p> <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • KUP Color Doppler (US only)
v7.1	<p>Torso/Torso-One Features:</p> <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • Trio • AI-assisted EF Workflow • PW Doppler • CW Doppler • TDI • AI FAST* • Cardiac Calcs • DICOM SR • KUP TDI mode (US only) <p>Lexsa Features:</p> <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • KUP Color Doppler (US only)
v7.2	<p>Torso/Torso-One Features:</p> <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • Trio† • AI-assisted EF Workflow • PW Doppler • CW Doppler • TDI

	<ul style="list-style-type: none"> • AI FAST*† • Cardiac Calcs • DICOM SR • KUP TDI mode (US only) • DICOMDIR support <p>Lexsa Features:</p> <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • KUP Lexsa Color Doppler
--	---

* Please note that although AI FAST is included in the software update, new customers will need to purchase the AI FAST license for each probe to activate the AI FAST feature.

† For the EU, Kosmos Trio and AI FAST are for educational use only.

Kosmos on Android: Purchasable Software Features

Software Version	Features
v1.0	<p>Torso-One Features:</p> <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • Trio • AI-assisted EF Workflow
v1.0.3 (Canada)	<p>Torso-One Features:</p> <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • AI-assisted EF Workflow
v2.0	<p>Torso-One Features:</p> <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • Trio • AI-assisted EF Workflow

	Lexsa Features: <ul style="list-style-type: none"> • B-mode • M-mode
v3.0	Additional Release Features: <ul style="list-style-type: none"> • Expansion of compatible tablets Torso-One Features: <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • Trio • AI-assisted EF Workflow Lexsa Features: <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler
v4.0	Torso-One Features: <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • AI-assisted EF Workflow • PW Doppler* • CW Doppler* • TDI* • Trio* • AI FAST* Lexsa Features: <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler
v4.1	Integrated Software Features <ul style="list-style-type: none"> • Kosmos UP • Us2.ai • 19Labs Torso-One Features: <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • AI-assisted EF Workflow • PW Doppler* • CW Doppler* • TDI* • Trio*†

	<ul style="list-style-type: none"> • AI FAST*† Lexsa Features: <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler
v5.1	Integrated Software Features <ul style="list-style-type: none"> • Kosmos UP • Us2.ai • 19Labs Torso-One Features: <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • AI-assisted EF Workflow • PW Doppler* • CW Doppler* • TDI* • Trio* • AI FAST* • Cardiac Calcs (Manual) • DICOM SR Lexsa Features: <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • Color Power Doppler • Vascular Calcs (Manual) • PW Doppler*

*** Please note that CW Doppler, PW Doppler, TDI, Trio, and AI FAST are licensed features. For new Android customers, licensed features can be purchased as yearly subscriptions or perpetual licenses.**

† For the EU, Kosmos Trio and AI FAST are for educational use only.

Kosmos on iOS: Purchasable Software Features

Software Version	Features
v1.1	<p>Torso-One Features:</p> <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler <p>Lexsa Features:</p> <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler
v2.1	<p>Torso-One Features:</p> <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • AI-assisted EF Workflow • PW Doppler* • CW Doppler* • TDI* • Trio* • AI FAST* • Cardiac Calcs (Manual) • DICOM SR <p>Lexsa Features:</p> <ul style="list-style-type: none"> • B-mode • M-mode • Color Doppler • Color Power Doppler • Vascular Calcs (Manual) • PW Doppler*

* Please note that CW Doppler, PW Doppler, TDI, Trio, and AI FAST are licensed features. For new Kosmos on iOS customers, licensed features can be purchased as yearly subscriptions or perpetual licenses.