

## **Kosmos on Android System Software v1.0.3**

Date: July 31<sup>st</sup>, 2023

Distribution: Internal (sales and clinical educators) and External (customers)

### **Kosmos on Android Software v1.0.3 Release Features**

The v1.0.3 system software release provides the following features for Kosmos on Android:

- Release of Kosmos on Android with Torso-One functionality
  - B-mode
  - M-mode
  - Color Doppler
  - EF Workflow

**Please note that Kosmos Trio is not available for use in Canada.**

### **Kosmos on Android v1.0.3 Software Availability**

The v1.0.3 system software is available for customers in Canada.

Visit the Google Play Store to download and install the latest update for the EchoNous Kosmos Ultrasound app. If you have any questions, contact EchoNous Customer Service by telephone at (844) 854-0800 or via email at [service@echonous.com](mailto:service@echonous.com).