



LEXSA™





A TRUE HYBRID. CART-BASED PERFORMANCE; HAND-CARRIED DEVICE.

EchoNous™ has developed a diagnostic platform aimed at shattering the paradigm of the use of ultraportable ultrasound devices in medicine, removing the diagnostic limitations usually associated with hand carried or quick look products. The Kosmos platform has advanced diagnostic and quantification capabilities not available on any other POCUS tool in its class, and uses AI and deep learning algorithms to facilitate training and automate complex calculations.

Continuing with changing perceptions, Lexsa is a 64/128 channel linear probe that uses the power of the Kosmos engine to produce outstanding image clarity for nerve, vascular, MSK, and lung applications.

Built for real-life patient situations

Patients come in all shapes and sizes, presenting with various pathologies. Kosmos was built with this in mind and tested to ensure high-confidence imaging on the most challenging patients. The proof is in the imaging.




One Probe, Multiple Configurations



Standalone Probe

The Lexsa probe is specially optimized to provide exceptional linear array imaging for MSK, Vascular, Nerve, and Lung applications.



Kosmos Mobile

Tailored for on-the-go clinicians, Kosmos Mobile includes an Apple iPad Air or Pro®, the Kosmos Essentials Carrying Case, Lexsa, Torso-One, or both probes, advanced Doppler features, and our AI capabilities.



Kosmos Plus

Fully-featured POCUS configuration with everything included – cutting-edge AI, advanced Doppler features, Lexsa, Torso-One, or both probes, a 13-inch Apple iPad®, and a medical-grade stand.

LEXSA SPECS

- Exam Types: Vascular, Nerve, MSK, Lung
- Modes: B-Mode, M-Mode, PW, Color, CPD
- 2D Measurements: Caliper, Distance
- Scan Depth: 1-10 cm
- 2D Frequency: 3-10.5 MHz
- Center Frequency: 7.5 MHz
- Footprint: 38mm x 8mm
- Mass: 280g
- Cable Length: 1.5 m



Durable

Lexsa is IP67 rated and tested to withstand the rigors of any point-of-care environment

5

5-Year Warranty

Industry-leading, no fine-print, 5-year warranty




Made in USA

Manufactured with excellence in Washington State

SEE KOSMOS IN ACTION

 (844) 854-0800
www.echonous.com

 **EchoNous Inc.**
8310 154th Ave NE
Building B, Suite 200
Redmond, WA 98052, U.S.A.

*The Kosmos Trio tool, AI FAST, and AI-Assisted EF workflow are released under the 'Enforcement Policy for Imaging Systems During the Coronavirus Disease 2019 (COVID-19) Public Health Emergency, Guidance for Industry and Food and Drug Administration Staff, April 2020' policy and have not been cleared by the FDA. These features are intended to be used by qualified healthcare professionals or under the supervision or in-person guidance of a trained or licensed healthcare professional.

© 2015 to 2024 EchoNous, Inc., or its affiliates. All rights reserved. EchoNous, echonous.com, and the EchoNous logo are trademarks of EchoNous, Inc. or its affiliates. EchoNous, Inc. 8310 154th Ave. NE, Building B, Suite 200, Redmond, WA 98052. D009673_Rev L