



Kosmos Bladder

Kosmos Bladder


NEXT GENERATION BLADDER VOLUME SOLUTION

Powered by the advanced AI technology of Kosmos, the Kosmos Bladder solution sets a new standard for accuracy, reliability, and user-friendly design—all while remaining affordable for any budget. Use it as a dedicated bladder scanner or expand its capabilities for an all-in-one POCUS solution.



Schedule a demo of Kosmos Bladder

 (844) 854-0800
www.echonous.com

 EchoNous Inc.
8310 154th Ave NE
Building B, Suite 200
Redmond, WA 98052, U.S.A.

Scan QR code



Kosmos Torso-One Transducer



Kosmos Mobile



All-in-one POCUS Solution with anti-theft lock

Fully detachable for ultimate flexibility



Kosmos Bladder represents a groundbreaking ultrasound solution that sets a new standard for bladder volume measurement.



Accurate



Reliable



All-in-one

Clinically Validated Accuracy

Kosmos Bladder AI can detect volumes from 4mL to over 1400mL and sets a new benchmark for accuracy, delivering precise measurements with ± 3 mL for volumes under 100 mL and $\pm 3\%$ for volumes over 100 mL.

Reliability You Can Count On

Kosmos is made in the USA with quality you can trust. Backed by a robust 5-year warranty, our military-grade transducers can withstand 1m drops and extreme environmental conditions.

More Than Just A Bladder Scanner

Expand Kosmos Bladder into a comprehensive AI-driven Point-of-Care Ultrasound solution. Offering bladder, vascular access, and multiple POCUS applications, Kosmos Bladder delivers an easy-to-use, versatile ultrasound device designed to empower healthcare professionals with accurate & reliable diagnostic capabilities.

*Not all features are available in all markets, please check with your local representative for availability in your region.