

# KOSMOS™ PLUS

## THE MOST POWERFUL, VERSATILE, AFFORDABLE POCUS SOLUTION FOR CLINICIANS EVERYWHERE.

Experience Kosmos Plus and unleash advanced POCUS capabilities at a fraction of the cost. With everything included — cutting-edge AI, advanced Doppler features, transducers, a 13-inch Apple iPad® and medical-grade stand — you'll have what you need to get started out of the box, and ensure your equipment stays together when you need it. Enhance patient care, optimize your workflow, and elevate your clinical information.



### Configuration

Probes, 13 inch iPad® Pro, medical-grade stand, AI software and advanced features all included



### Durable

Probes designed for the point of care: survive 1 m drop, IP67 rated, and approved for all levels of cleaning and disinfection



### 5-Year Warranty

Industry-leading, no-fine-print 5-year warranty for all transducers



### Made in USA

Manufactured with excellence in Washington State



# Kosmos Plus is designed for you:



Cardiology



Critical Care



Emergency  
Medicine



Vascular  
Access







Primary  
Care



Nephrology

## AI + Doppler Features

Our intelligent POCUS solution on iOS provides AI-assisted workflows (AI FAST, AI-Assisted EF Workflow, Auto Doppler, and Auto Preset) and advanced Doppler imaging (PW, TDI, CW, Color Doppler, Color Power Doppler), to give you the clinical insights you need for every care scenario.



For the first time, you'll have the capabilities you need for the most complicated cases balanced with the versatility, portability, and affordability of a handheld system.

### Stand

- Weight: 24.5 kg
- Dimensions: 112 x 51 x 45 cm

### Torso-One™


- Weight: 275 g
- Footprint: 20.7 mm x 16 mm
- 2D Frequency: 1.5-4.5 MHz
- Center Frequency: 3 MHz

### Lexsa™

- Weight: 280 g
- Footprint: 38 mm x 8 mm
- 2D Frequency: 3-10.5 MHz
- Center Frequency: 7.5 MHz

## SEE KOSMOS IN ACTION

 (844) 854-0800  
[www.echonous.com](http://www.echonous.com)

 **EchoNous Inc.**  
8310 154th Ave NE  
Building B, Suite 200  
Redmond, WA 98052, U.S.A.

%Not all features are available in all markets, please check with your local representative for availability in your region.

© 2015 to 2025 EchoNous, Inc., or its affiliates. All rights reserved. EchoNous, echonous.com, and the EchoNous logo are trademarks of EchoNous, Inc. or its affiliates. EchoNous, Inc. 8310 154th Ave. NE, Building B, Suite 200, Redmond, WA 98052. D011210\_REV F