



KOSMOS LINK



UNINTERRUPTED POWER FOR UNINTERRUPTED CARE

Kosmos Link enables continuous operation with charge-while-scanning capability, multi-probe connectivity, and enhanced power delivery for advanced imaging.

Offering hours of scanning on battery power alone—a previously unheard of amount of time in the portable ultrasound world—Kosmos Link adapts to your clinical environment, whether deployed as a tabletop or cart configuration.



Charge while
Scanning



Multi-Probe
connectivity



Extended Battery
Life



Cart to Handheld
Portability

Kosmos Link Works Where You Do




On a Stand
Anti-theft lock and
quick release



On a Table Top or Bed
Built-in kickstand for a table
top or the bed



On an IV Pole
Built-in design for ability to
hang on an IV Pole



Technical Specifications



Dimensions:

292.16 mm x 224.5 mm x 34.35 mm



Weight:

1050 g



Battery Life:

A fully charged Kosmos Link will provide 3-8 hours of uninterrupted scanning (performance varies based on scanning modes used).



Battery Capability:

Ability to charge while scanning



Battery Charge Time:

0-90% in ~2 hours



Multiple Ports:

3 ports to support 3 probes



Compatible Tablets:

Please visit: <https://echonous.com/product/device-compatibility/> for tablet compatibility



"Kosmos Link has had an immediate impact on our echocardiogram workflow. Allowing us to quickly and easily shift between different study acquisitions without having to unplug and plug each probe... The built-in battery and ability to keep the device charging while we're using it means we can do more studies in a day."

Henry L. Burkholder, IV MD



 (844) 854-0800
 www.echonous.com

 **EchoNous Inc.**
8310 154th Ave NE
Building B, Suite 200
Redmond, WA 98052, U.S.A.

*Not all features are available in all markets, please check with your local representative for availability in your region.

EchoNous, Inc., or its affiliates. All rights reserved. EchoNous, echonous.com, and the EchoNous logo are trademarks of EchoNous, Inc. or its affiliates. D012184_Rev A