




KOSMOS LINK

UNINTERRUPTED POWER FOR UNINTERRUPTED CARE

Kosmos Link enables continuous operation with charge-while-scanning capability, multi-probe connectivity, and enhanced power delivery for advanced imaging.

Offering hours of scanning on battery power alone—a previously unheard of amount of time in the portable ultrasound world—Kosmos Link adapts to your clinical environment, whether deployed as a tabletop or cart configuration.



Charge while
Scanning



Multi-Probe
connectivity



Industry-Leading
Battery Life



Cart-to-Handheld
Portability

Kosmos Link Works Where You Do



On a Stand

Anti-theft lock and
quick release




On a Table Top or Bed

Built-in kickstand for a table
top or the bed



On an IV Pole

Built-in design for ability to
hang on an IV Pole



Technical Specifications



Dimensions:

8.85 x 11.6 x 1.125 in



Weight:

950 g

Battery Life:



A fully charged Kosmos Link + iPad Pro® will provide up to 2.8 hours of continuous B-mode scanning (performance varies based on iPad model and scanning modes used)*



Battery Capability:

Ability to charge while scanning



Battery Charge Time:

0-90% in ~2 hours



Multiple Ports:

3 ports to support 3 probes



Compatible Tablets:

Please visit: <https://echonous.com/product/device-compatibility/> for tablet compatibility



"Kosmos Link has had an immediate impact on our echocardiogram workflow. Allowing us to quickly and easily shift between different study acquisitions without having to unplug and plug each probe... The built-in battery and ability to keep the device charging while we're using it means we can do more studies in a day."

Henry L. Burkholder, IV MD



(844) 854-0800
www.echonous.com

EchoNous Inc.
 8310 154th Ave NE
 Building B, Suite 200
 Redmond, WA 98052, U.S.A.

*Not all features are available in all markets, please check with your local representative for availability in your region.

EchoNous, Inc., or its affiliates. All rights reserved. EchoNous, echonous.com, and the EchoNous logo are trademarks of EchoNous, Inc. or its affiliates. D012184_Rev C