



LEXSA®




A TRUE HYBRID. CART-BASED PERFORMANCE; HAND-CARRIED DEVICE.

EchoNous™ has developed a diagnostic platform aimed at shattering the paradigm of the use of ultraportable ultrasound devices in medicine, removing the diagnostic limitations usually associated with hand carried or quick look products. The Kosmos platform has advanced diagnostic and quantification capabilities not available on any other POCUS tool in its class, and uses AI and deep learning algorithms to facilitate training and automate complex calculations.

Continuing with changing perceptions, Lexsa is a 64/128 channel linear probe that uses the power of the Kosmos engine to produce outstanding image clarity for MSK, Nerve, Vascular, Lungs, & VAF applications.

Built for real-life patient situations

Patients come in all shapes and sizes, presenting with various pathologies. Kosmos was built with this in mind and tested to ensure high-confidence imaging on the most challenging patients. The proof is in the imaging.



One Probe, Multiple Configurations



Standalone Probe

The Lexsa transducer is specially optimized to provide exceptional linear array imaging for MSK, Nerve, Vascular, Lung, and VAF applications.



Kosmos Mobile

Tailored for on-the-go clinicians, Kosmos Mobile includes an Apple iPad®, Kosmos Essentials Carrying Case, Lexsa, Torso-One, or both transducers, Doppler features, and AI capabilities.



Kosmos Plus

Fully-featured POCUS configuration with everything included – AI, Doppler features, Lexsa, Torso-One, or both transducers, choice of a Apple iPad® or Kosmos Tablet G1, and a medical-grade stand.

LEXSA SPECS

- Exam Types: MSK, Nerve, Vascular, Lung, & VAF
- Modes: B-Mode, M-Mode, PW, Color, CPD
- 2D Measurements: Caliper, Distance
- Scan Depth: 1-10 cm
- 2D Frequency: 3-10.5 MHz
- Center Frequency: 7.5 MHz
- Footprint: 38mm x 8mm
- Mass: 280g
- Cable Length: 1.5 m



Durable

Lexsa is IPX7 rated and tested to withstand the rigors of any point-of-care environment



5-Year Warranty

Industry-leading, no fine-print, 5-year warranty



Made in USA

Manufactured with excellence in Washington State

SEE KOSMOS IN ACTION



(844) 854-0800
www.echonous.com



EchoNous Inc.
8310 154th Ave NE
Building B, Suite 200
Redmond, WA 98052, U.S.A.

*Not all features are available in all markets, please check with your local representative for availability in your region.

© EchoNous, Inc., or its affiliates. All rights reserved. EchoNous, echonous.com, and the EchoNous logo are trademarks of EchoNous, Inc. or its affiliates. EchoNous, Inc. 8310 154th Ave. NE, Building B, Suite 200, Redmond, WA 98052. D009673_Rev T