




KOSMOS® PLUS

THE MOST POWERFUL, VERSATILE, AFFORDABLE POCUS SOLUTION FOR CLINICIANS EVERYWHERE.

Experience Kosmos Plus and unleash advanced POCUS capabilities at a fraction of the cost. With everything included — cutting-edge AI, advanced Doppler features, flexible transducers, Apple iPad® or Kosmos Tablet G1, and medical-grade stand — you'll have what you need to get started out of the box, and ensure your equipment stays together when you need it. Enhance patient care, optimize your workflow, and elevate your clinical information.



Configuration

Flexible transducers, Apple iPad® or Kosmos Tablet G1, medical-grade stand, AI software and advanced features all included



Durable

Kosmos transducers are drop tested from 1m, IPX7 rated, and have robust chemical compatibility



5-Year Warranty

Industry-leading, no-fine-print 5-year warranty for all transducers




Made in USA

Manufactured with excellence in Washington State



Why Choose Kosmos Plus?

Kosmos Plus gives you the performance of a fully-featured POCUS cart at a fraction of the cost. With advanced AI (AI FAST, Auto EF, Auto Doppler, and Auto Preset), Doppler capabilities (PW, TDI, CW, Color Doppler, Color Power Doppler), and flexible transducer options — including Torso-One for whole torso imaging and Lexsa for high-resolution MSK, Nerve, Vascular, Lung, and Vascular Auto Flow applications — you're equipped to handle a wide range of clinical needs with confidence.



Kosmos Stand

- Weight: 24.5 kg
- Dimensions: 16 x 23 x 53in

Kosmos Link

- Weight: 950 g
- Dimensions: 8.85 x 11.6 x 1.125in



Torso-One™

- Phased Array Probe
- Presets: Heart, Lung, Abdomen, OB, GYN, & Bladder
- Weight: 275 g
- Footprint: 20.7 mm x 16 mm
- 2D Frequency: 1.5-4.5 MHz
- Center Frequency: 3 MHz



Lexsa

- Linear Array Probe
- Presets: MSK, Nerve, Vascular, Lung & VAF
- Weight: 280 g
- Footprint: 38 mm x 8 mm
- 2D Frequency: 3-10.5 MHz
- Center Frequency: 7.5 MHz

SEE KOSMOS IN ACTION

☎ (844) 854-0800
www.echonous.com

📍 **EchoNous Inc.**
8310 154th Ave NE
Building B, Suite 200
Redmond, WA 98052, U.S.A.

Note: Not all features are available in all markets, please check with your local representative for availability in your region.

EchoNous, Inc., or its affiliates. All rights reserved. EchoNous, echonous.com, and the EchoNous logo are trademarks of EchoNous, Inc. or its affiliates. EchoNous, Inc. 8310 154th Ave. NE, Building B, Suite 200, Redmond, WA 98052. D011210_Rev K