



Kosmos Tablet G1

Built to Perform in the Most Demanding Environments

Custom-engineered with Getac, Kosmos Tablet G1 redefines rugged POCUS. By merging military-grade durability with a streamlined user experience, it provides a powerful, reliable imaging solution specifically built for the high-stake demands of diverse clinical environments.

The Kosmos Tablet G1 is designed to work as hard as you do, right out of the box.

- **Purpose-Built:** Ready to use out of the box with the Kosmos platform pre-installed
- **Rapid Startup:** Boots up in under 30 seconds – just connect to Wi-Fi to register your probe and start scanning
- **Enterprise-Ready:** Dramatically reduces IT setup time for easy deployment
- **The Result:** Less time troubleshooting technology, and more time delivering critical patient care



 **Touchscreen Display**
10-inch fully rugged interface

 **Volume Control**
Easily adjust system audio on the fly

 **Charging Indicator**
Confirm charging status via green-amber lighting at a glance

 **Managing Power**
Various power options for zero downtime

Technical Specifications



Software:

Pre-loaded with the Kosmos software



Dimensions:

275 (X) x 195 (H) x 17.9 (D) mm



Screen Size:

10.1 inches



Weight:

1020.58 grams



Display:

Touchscreen Display



Battery Life:

150 minutes continuous scanning
(240 minutes with Kosmos Link)

"Kosmos Tablet G1 is the platform our most demanding customers have been asking for. By partnering with Getac, we've been able to deliver a custom-engineered, hospital-grade Android solution that is as rugged as the environments it's deployed in and as intelligent as the Kosmos platform it runs. From the battlefield to the bedside, Kosmos Tablet G1 is ready to go from the moment it's unboxed."

- Tony Titus, CEO of EchoNous



Handstrap/Kickstand

Secure strap for one-handed mobility or hands-free kickstand



Kosmos AI Station

- Medical-Grade Stand
- Weight: 24.5 kg
- Dimensions: 16 x 23 x 53 in



Kosmos Link

- Charge-while-scanning capability
- Weight: 950 g
- Dimensions: 8.5 x 11.6 x 1.125 in

SEE KOSMOS IN ACTION

(844) 854-0800

www.echonous.com

EchoNous Inc.
8310 154th Ave NE
Building B, Suite 200
Redmond, WA 98052, U.S.A.



*Not all features are available in all markets, please check with your local representative for availability in your region.

EchoNous, Inc., or its affiliates. All rights reserved. EchoNous, echonous.com, and the EchoNous logo are trademarks of EchoNous, Inc. or its affiliates.
D012553_Rev B