

Your NXT Generation Transducer **Torso NXT**

Torso NXT is our next-generation ultrasound transducer, built to extend point-of-care imaging across clinical settings. Bringing cart-level diagnostic capability to an ultraportable probe with a 360° ergonomic design, Torso NXT combines comfort, dependability, and clinical precision on the Kosmos platform.

WHOLE-TORSO IMAGING

Heart, lung, abdomen, OB/GYN and bladder—from a single probe.



World-class Imaging

Experience the world-class image quality and capabilities of a premium ultrasound cart in the palm of your hand, driven by the Torso NXT's proprietary acoustic tuning and PZT hardware.

Comfort, redefined

A contoured 360° grip and a 50% reduction in volume may help reduce scanning fatigue during high-volume shifts.

Minimal downtime

A user-replaceable cable with a reinforced threaded collar. Swap it out on-site for maximum uptime.



IP67 rated

Tested to withstand the rigors of point-of-care environment.



Up to 5-year warranty*

2- to 5-year coverage, with optional extensions available.




Made in the USA

Manufactured with excellence in Washington State.

Cart-Level Imaging **in Your Hand.**

Kosmos delivers image quality and performance comparable to cart-based systems, with the portability and affordability of a handheld.



One probe. Multiple configurations.



Standalone Probe

The Torso NXT probe on its own—cart-level, whole-torso imaging for heart, lung, abdomen, OB/GYN and bladder.



Kosmos Mobile

For on-the-go clinicians: Apple iPad, Kosmos Essentials carrying case, your choice of Torso NXT or Lexsa transducer, Doppler and AI.



Kosmos Plus

Fully featured POCUS: AI, Doppler, your choice of probe, Apple iPad or Kosmos Tablet G1, and a medical-grade stand.

TECHNICAL SPECIFICATIONS

TRANSDUCER

PZT phased array

SCAN DEPTH

4–30 cm

FREQUENCY RANGE

1.5–4.5 MHz

CENTER FREQUENCY

3 MHz

FOOTPRINT

20.7 × 16 mm

WEIGHT

284 g

SIZE L × W × D

163 × 52 × 35 mm



BOOK A DEMO
(844) 854-0800

www.echonous.com

EchoNous, Inc. · 8310 154th Ave NE, Building B, Suite 200, Redmond, WA 98052, U.S.A.

Not all features are available in all markets. Check with your local representative for availability in your region.

*Warranty terms and optional extensions are region-specific; consult your local representative for details.

© EchoNous, Inc., or its affiliates. All rights reserved. EchoNous, echonous.com, and the EchoNous logo are trademarks of EchoNous, Inc. or its affiliates.

D012886_Rev A